

# 941

## SENSORS AND ACTUATORS

### NURSE CALL & PRESENCE BUTTON



Code : **9410054**

Model : **LLC-142**

#### Description

Room peripheral with two buttons: a nurse call button with a red LED and a nurse presence button with a green LED. This peripheral allows patients to call the nurse with the red button. Once initiated the call, the red LED helps to locate the origin of the call easily. The green button is used to either indicate that the nurse is present or to cancel the call. Wall or surface mount.

9410054

LLC-142

#### Especificaciones

Botón de llamada a enfermera	Color botón: Rojo (RAL 3020, rojo tráfico) Color indicador led: Rojo Icono: enfermera Color icono: blanco (RAL 9003, blanco señales)
Botón de presencia de enfermera	Color botón: verde (RAL 6018, verde amarillento) Color indicador led: verde Icono: persona 1 Color icono: blanco (RAL 9003, blanco señales)
Alimentación	24 V
Consumo (reposo)	10 mA
Consumo (máximo)	40 mA
Dimensiones	80 x 80 x 30 mm
Peso	80 g
Material	ABS / Policarbonato
Color (embellecedor)	Blanco (RAL 9001, blanco crema)
Temperatura de funcionamiento	5°C a 55°C
Protección	IP 42
Montaje	Empotrado (con ACC-100 / ACC-101) Superficial (con ACC-110)
Normativa (Seguridad)	EN 60950-1:2006
Normativa (Emisiones electromagnéticas)	EN 61000-3-3

